



Model #: ATN-R06N



Overview

An optional conversion kit for adapting an Attenda bed for use as a restraint bed. Conversion kit consists of six heavy-duty restraint rings that mount through hole locations in bed and through-bolt into floor. Restraint rings then serve as connector points for user-supplied restraint devices. Norix does not provide the restraint devices.

Material:

1/4" thick steel.

Weight:

6 Rings: 1.1 lbs.

Assembly:

None required.

Installation:

It is intended that the restraint ring is installed with the same anchors used to secure the bed. Anchors not supplied by Norix.

Intended Use Disclaimer:

These goods are sold by Norix for institutional use only and not as consumer products. Except for the warranty of title, NO WARRANTY OF MERCHANTABILITY, FITNESS OR OTHER WARRANTY (WHETHER EXPRESSED, IMPLIED OR STATUTORY) IS MADE BY NORIX, except that Norix warrants the goods to be free from defects in materials and workmanship in normal use and service. NORIX IS NOT LIABLE FOR CONSEQUENTIAL DAMAGES.

THESE GOODS ARE NOT WARRANTED AS SUITABLE FOR ANY PARTICULAR PURPOSE PARTICULAR TO BUYER. THE SUITABILITY OF GOODS FOR ANY PURPOSE PARTICULAR TO BUYER IS FOR BUYER, IN BUYER'S SOLE JUDGEMENT, TO DETERMINE. NORIX ASSUMES NO RESPONSIBILITY FOR THE SELECTION OR FURNISHING OF GOODS SUITABLE TO THE INDIVIDUAL NEEDS AND PURPOSES OF ANY PARTICULAR BUYER.

Warranty:

5 year limited replacement warranty.

Maintenance:

Easily cleaned with a damp cloth and a mild soap solution.

Filing System:

12600 Institutional Furniture
11197 Correctional Furniture